


Medtronic

# Cordless freedom to do more,<sup>1,†,‡, ΩΩ</sup> your way.

The next advancement<sup>†</sup> in cordless ultrasonic dissection is here – now designed to coagulate and dissect vasculature up to and including 7 mm in diameter.<sup>1,‡</sup>



## A cordless design facilitates:

- More natural movement<sup>2,§</sup>
- Easier instrument exchanges<sup>2,Ω</sup>
- Improved use of OR space<sup>3,††</sup>

## A tapered, curved jaw provides:

- Access in tight spaces<sup>2,††</sup>
- Visualization of target structures<sup>2,††</sup>
- Precise energy dissection<sup>1,‡</sup>

## A single intuitive<sup>1,‡</sup> dual mode energy activation button enables:

- Procedural versatility<sup>1,§§</sup>
- The surgeon's eyes to stay on the surgical field<sup>1,ΩΩ</sup>
- Improved surgical efficiency<sup>1,‡</sup>

# The cordless freedom you expect with the advancements<sup>†</sup> you desire.



The Sonicision™ 7 curved jaw device remains the first and only cordless ultrasonic dissector. Now with the ability to coagulate and dissect vessels up to and including 7 mm in diameter.<sup>1,‡</sup>

## Cordless power in the palm of your hand

A reusable generator enables:

- Use for **150 procedures** before replacement<sup>4</sup>
- **Streamlined<sup>†</sup> cleaning and sterilization** because it's autowashable and autoclavable

A reusable nonsterile battery pack provides:

- **88 percent greater capacity<sup>5,†</sup>** to instill confidence of consistent power throughout long procedures
- **200 procedural uses<sup>6,7</sup>** – twice the number of earlier devices<sup>†</sup>
- No-hassle power management via **a simple, intuitive charger<sup>3,†††</sup>**



## Ordering information

|                |  |        |
|----------------|--|--------|
| <b>SCD7A13</b> | Sonicision™ 7 curved jaw cordless ultrasonic dissector – 13 cm | 6/case |
| <b>SCD7A26</b> | Sonicision™ 7 curved jaw cordless ultrasonic dissector – 26 cm | 6/case |
| <b>SCD7A39</b> | Sonicision™ 7 curved jaw cordless ultrasonic dissector – 39 cm | 6/case |
| <b>SCD7A48</b> | Sonicision™ 7 curved jaw cordless ultrasonic dissector – 48 cm | 6/case |
| <b>SCBA</b>    | Sonicision™ reusable battery pack                              | 1/each |
| <b>SCG7AB</b>  | Sonicision™ 7 reusable generator                               | 1/each |
| <b>SCBIGA</b>  | Sonicision™ reusable battery insertion guide                   | 2/pack |
| <b>CBCA</b>    | Sonicision™ battery charger                                    | 1/each |
| <b>SCSTA</b>   | Sonicision™ reusable sterilization tray                        | 1/each |

<sup>†</sup>Compared to earlier generations of the Sonicision™ curved jaw device. <sup>‡</sup>25 out of 27 surgeons surveyed after use agreed. <sup>§</sup>30 out of 33 surgeons surveyed after use agreed. <sup>¶</sup>32 out of 33 surgeons surveyed after use agreed. <sup>††</sup>29 out of 29 nurses surveyed after use agreed. <sup>†††</sup>33 out of 33 surgeons surveyed after use agreed. <sup>§§</sup>27 out of 27 surgeons surveyed after use agreed. <sup>¶¶</sup>26 out of 27 surgeons surveyed after use agreed. <sup>††††</sup>28 out of 28 nurses surveyed after use agreed. <sup>†††††</sup>29 out of 33 surgeons surveyed after use agreed. **1.** Based on internal test report #RE00329878 rev A, Marketing evaluation of surgeon experience using the Sonicision™ 7 curved jaw cordless ultrasonic dissector. Apr. 14-15 and 20-22, 2021. **2.** Based on internal test report #R0042752 rev A, Marketing evaluation of surgeon experience using the Sonicision™ curved jaw cordless ultrasonic dissector. Feb. 13-17, Feb. 28, March 3, May 16, 2017. **3.** Based on internal test report #R0049393 rev A, Operating room staff marketing evaluation of the Sonicision™ curved jaw cordless ultrasonic dissector. Dec. 1, 6-7, 2016. **4.** Based on internal test report #RE00108733, rev A, Reusable life test verification report (SCGAA-TA9). August 24, 2017. **5.** Based on internal test memo #RE00153542, rev B, Battery capacity comparison between SCB and SCBA. May 17, 2018. **6.** Based on internal test report #R0042707, rev A, Reusable life test verification report for the SCBA reusable battery pack. April 26, 2017. **7.** Based on internal test report #RE00113195, rev A, Reusable life test verification report for the SCBA reusable battery pack – Oxivir. September 13, 2017.

Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2022 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. <sup>™\*</sup> Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. EMEA-SE-2200095 - 22-weu-sonicision-7-info-sheet-7269379

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)

Bring greater efficiency<sup>1,‡</sup> and freedom of movement<sup>2,†††</sup> to your OR.  
Call your Medtronic sales rep today or visit [medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)

**Medtronic**