



Surgipro™ suture

Surgipro™ II suture

Reference Guide



Surgipro™ and Surgipro™ II suture is a brand you have trusted for many years with several product innovations.

Timeline

- 1991** Surgipro™ suture launched
- 1994** Surgalloy™ premium needle alloy introduced
- 1996** Lubricant added to help minimize suture breakage
- 2001**
 - Surgipro™ II suture launched with improved handling characteristics
 - NuCoat™ coating added to improve needle with multipass penetration
- 2007**
 - Surgipro™ II suture resin enhanced to improve knot run down characteristics
 - NuPack™ designed to reduce memory in suture
- 2021**
 - Sharper, stronger Cardiopoint™ needles
 - Strand diameter optimization
 - NuPack™ retainer improvement

Frequent Uses

- Soft tissue approximation
- Ophthalmic
- Cardiovascular
- Neural tissue

Surgipro™ II Suture Top Codes

Code	Color	Size USP	Length		Needles
			cm	Inches	
VP-709-X	Blue	6-0*	75	30	CV-1
VP-557-X	Blue	4-0*	90	36	CV-23
VP-703-X	Blue	6-0*	60	24	CV-1
SP-5698-G	Blue	5-0	45	18	P-13
VP-702-X	Blue	7-0*	60	24	CV-1

Product Comparison

Product	Surgipro™ Suture ¹ Surgipro™ II Suture ¹	Prolene™ Suture ²
Structure	Monofilament	Monofilament
Suture Type	Non-Absorbable	Non-Absorbable
Composition	Polypropylene	Polypropylene
Coating	Uncoated	Uncoated
Color	Blue, Clear	Blue, Clear
Tensile Strength	Permanent	Permanent
Absorption Profile	Permanent	Permanent
Sizes	Surgipro™ suture: Sizes 10-0, 6-0, 2-0 through 2 Surgipro™ II suture: Sizes 8-0 to 3-0	10-0 through 2
Indications	Surgipro™ and Surgipro™ II polypropylene suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neural tissue.	Prolene™ polypropylene suture is indicated for use in general soft tissue approximating and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.
Contraindications	None known	None known
Box Quantities	1, 2 or 3 dozen	1, 2 or 3 dozen

*Double-arm

1. Surgipro™ II Instructions for Use.

2. Prolene™ Instructions for Use.

© 2021 Medtronic. All rights reserved. EMEA-WC-2100048 - 21-weu-surgipro-product-guide-5502633

Medtronic
Further, Together